

## LABORATORY EQUIPMENT Laboratory Sieve

### Laboratory Sieve Type LN for Wet Screening

The apparatus for wet analyses consists of laboratory apparatus Type L with Transformer 220/24 V; two stands; inlet funnel with tap 10 mm, of glass, capacity 500 cm<sup>3</sup>; collecting pan; holder for five screen trays 100 mm Ø open top and bottom.

This apparatus is specially suitable for wet analyses of any kind, particularly of sludges, etc. Absolutely accurate results can be obtained by washing with water or other fluids.



### Laboratory Sieve Type E as Installation Unit

The "Schallfix" laboratory sieve can also be supplied as an installation unit. In that construction the lower connection is in the form of a flange and the tappet end is in the form of a collet.

The E-apparatus can serve the most varied purposes, firmly fixed to tables, vessels, walls, etc. screens, funnels, etc., can be put into the collet. The E-apparatus is supplied for direct connection to 220 V.